

QuickSplit™ Multiport Fixed Flow Splitters

Divides the incoming flow stream into 3 or 4 channels. The *QuickSplit* Multiport Flow Splitter is ideal for applications that use multiple detectors and/or a fraction collector. Split ratios can be configured to be identical in each channel or custom configured to accommodate specific flow rates at each channel. The low internal dead volume prevents excessive dispersion and the replaceable inlet filter insures robust operation. Available in models for analytical and semi-preparative inlet flow rates. **Please contact ASI for assistance in configuring your *QuickSplit* Multiport Fixed Splitter prior to ordering. Please specify inlet flow rate and desired outlet flow rate in each channel when you order.**



***QuickSplit* Multiport Fixed Flow Splitter**