

Chemical Analyses of Tungsten Powder

Part#: 160-0012

1.2 μm Tungsten Particles, Pkg 50 gram

QUANTITATIVE ANALYSES					
Element	RT04000 Limits	Analytical Result	Element	RT04000 Limits	Analytical Result
ppm Al	10 ppm	4 ppm	FSSS μm	0.95-1.25	1.20 μm
Ca	50 ppm	<1 ppm	B.D. g/in ³	20-45	28.2 g/in ³
Cr	25 ppm	7 ppm	LOR	1200 ppm	1000 ppm
Cu	10 ppm	<1 ppm	Mo ppm	500 ppm	<8 ppm
Fe	60 ppm	14 ppm			
K	150 ppm	<10 ppm			
Mg	50 ppm	<1 ppm			
Mn	10 ppm	2 ppm			
Ni	100 ppm	7 ppm			
Si	50 ppm	2 ppm			
Sn	20 ppm	<1 ppm			
Na	100 ppm	5 ppm			

The above information covers the physical and chemical analyses of the Tungsten Powder shipped to you.